

Wearing Instructions for Shoulder Saddle Sling



To Wear Sling:

1. Place leather shoulder portion on top of affected shoulder.
2. Attach thick cotton strap across body and under opposite arm, hooking through latch.
3. Place leather forearm portion (with 2 velcro straps) around mid – forearm and tighten straps.
4. To support shoulder joint, pull velcro straps (that attach forearm leather to shoulder leather) tight enough to support arm against gravity.
5. To take off, simply un-velcro forearm portion and unlatch opposite arm strap.
6. Remove from shoulder

*Wear per your doctor's instructions.

*Wear sling on top of shirt to prevent skin breakdown.

*Should you have any questions or concerns, please consult your orthotist from Infinite Technologies O & P at contact numbers below.

